



Materials & Sundries

IMPRESSION MATERIALS

FILLING MATERIALS

CEMENTS AND LINERS

HYGIENE, ENDO AND POLISHING

INSTRUMENTS

www.sswwhite.com



SS White Manufacturing

Stephenson Drive, Gloucester, GL2 2HA, United Kingdom.

T: +44 (0) 1452 729 751 F: +44 (0) 1452 307 187

sales@sswhite.com www.sswwhite.com



Welcome to SS White

Since 1844 SS White have manufactured and distributed a range of quality professional dental products and are renowned the world round for their superior quality. SS White materials tradition continues today using the support of a network of specialist distributors to serve the dental professional all around the world and with this our latest catalogue showcasing.

SS White Manufacturing UK offers a particularly popular range of dental materials, including impression materials, waxes, cements, lining materials, cross infection control products and surgery sundries. The extensive SS White range offers reliable value for money materials for day to day dental practice.

The superiority of all SS White materials is assured by its robust quality systems as a part of our ISO 13485 accredited Quality Management System. In addition all products carry CE and FDA product registration.

Please also visit our partner Prima Dental Group for all your rotary and instrument requirements.



Impression Materials



Impression Paste

Zinc oxide and eugenol. The original Zinc Oxide and Eugenol Impression material. Unsurpassed for full impressions. Undergoes no dimensional change.


Code	Quantity
IP	Box of 2 tubes and mixing pad



Casting Wax Red No 6

Deep red wax for constructing patterns for casting saddles, bars, clasps, one-piece casting of dentures. Carves easily without chipping.

Code	Quantity
WX645	45 gm tin



Tenax Wax

High Beeswax Denture Wax. Medium soft, sheet size (143 mm x 75 mm) unique for toughness and pliability

Code	Quantity
WXTEN	227 gm (1/2 lb)
WXTEN5	2.27 kg (5 lb)



SS White (Willis) Bite Gauge

A simple measuring instrument as an aid in the re-establishment of correct facial contour during full or partial denture work. Slide B is fully adjustable and reversible along Shaft A and may be fixed in any position by means of its locking screw.

Filling Materials



SS White Advance 70

SS White is a high silver content, dispersed phased non gamma-2 amalgams alloy – with all material advanced of the latest amalgam technology. Optimum working properties offers the dentist as well as other characteristics of quality that assure excellent, corrosion resistant amalgam restorations with durable surface lustre.

Code	Quantity
ADVAN701	50 caps
ADVAN702	50 caps
ADVAN703	50 caps




Legend™ Silver

Glass Ionomer Posterior Filling Material

Code	Quantity
LS	15G

Legend™ Refills

Code	Quantity
Rapid A1, A2, A3.5, B2, B4, C3	12gm
Regular A2, A3.5, B2	12gm



New Advance 44.5% Silver

High copper, non-gamma, 2 amalgam, zinc free Powder

Product	Code	Quantity
Capsules	ADVAN1	50 caps
Capsules	ADVAN2	50 caps
Capsules	ADVAN3	50 caps



SS White Etching Gel

Fast action, excellent contrast with oral cavity, stays in place, will not "slump" and will not stain the tooth surface.

Code	Quantity
Green with 10 tips	ZZEGG 5ml
Green with 20 tips 12 ml pack	ZZEGG12 5 pack 12ml
Dispensing Tips	ZZEGN100 100



SS White New True Dentalloy

High copper, non-gamma, 2 amalgam, zinc free Powder

Code	Quantity
NTDA	50 gm
	250 gm

Cements and Liners



PCI Zinc Polycarboxylate Cement

Product	Code	Quantity
Powder	PCIP	60 gm
Liquid	PCIL	40 ml
Combination Pack	PCIPL	Each

Cements and Liners



ZCI (Zinc Cement Improved)

Zinc Oxyphosphate Cement. Ratio: 1.3g Powder to 0.5ml Liquid.
Colours: 11 Pure White, 12 Tooth yellow.

Product	Code	Quantity
Powders No 12	CZ...32	32 gm
Powders No 12	CZ...90	90 gm
Powders No 12	CZ...230	230 gm
Liquid	CZL15	15 ml
Liquid	CZL45	45 ml
Liquid	CZL125	125 ml

Hygiene, Endo and Polishing



Gobble Plus

Product	Code	Quantity
Evacuation system cleaner	GOBBLE2	2 ltr
Refill 10 litres	GOBBLE10	10 ltr

Gleam

Product	Code	Quantity
Ultrasonic cleaning solution	GLEAM2	2 ltr
Refill 10 Litres	GLEAM10	10 ltr

Hygiene, Endo and Polishing



Propylaxis Paste

Conventional pumice based peppermint flavoured excellent for stain removal and general purpose.

Code	Quantity
PP50	50 gm jar
PP200	200 gm jar
PP400	400 gm jar



Oraproph

Removes stubborn stains from tooth enamel with out damage - oil free.

Code	Quantity
PPOPH	200 gm jar



Mouthwash Tablets

Product	Code	Quantity
Pink. Conventional pink thymol flavoured	ZZMWP	Bottle of 1000
Green. Formulated with lime flavour to provide alternative	ZZMWG	Bottle of 1500